

## Protein Detection Using *Identity* Raman Plate Reader.

### INTRODUCTION



**Figure 1.** Digilab's Identity Raman Plate Reader. (Catalogue #RMI-78500-1 and #RMI-53200-1)

Biopolymer (proteins, nucleic acids, polyamines etc.) at high concentrations can be easily detected by Raman spectroscopy [1-4]. Unfortunately, most of natural samples have these substances present at significantly lower levels, making application of the conventional Raman spectroscopy towards natural samples disputable. Sensitivity of the biopolymer's detection can be significantly enhanced by application of Surface-Enhanced Raman Scattering (SERS) technique. This Application Note describes how to use Identity Raman Plate Reader in both conventional Raman, and SERS mode to successfully detect some proteins in solution at different concentrations, particularly at the levels below 1%.

### PROCEDURE


**Materials and Methods:** All measurements were performed with a Digilab Identity Raman plate reader (Figure 1). The system is configurable with either a 532 nm (Model #RMI-53200-1) or a 785 nm (Model #RMI-78500-1) laser. Raman scattering is collected in an 180° backscatter configuration by a spectrometer with a Peltier-cooled CCD array detector capable of  $<10\text{ cm}^{-1}$  spectral resolution. The Identity supports standard clear, flat bottom 96 and 384 well microtiter plates, as well as custom plate formats. The laser is focused through the bottom of the plate into the well for the analysis of liquids. Powders or coatings on the bottom of the well can be analyzed by placing the plate on a plate shim to raise the bottom of the well to the laser focal point. For SERS application, silver colloid was prepared according to Lee and Meisel [5] and was diluted 5-20 fold with corresponding buffer immediately before use. The colloid was mixed with the sample in 1:1 ratio inside the microplate well and measurement was taken within 1 min after the mixing and additionally 30 min later.

**Results and Discussion:** Conventional Raman spectroscopy has relatively low sensitivity (with lower detection limits being around 0.1%). Biopolymers in aqueous solutions can usually be detected only at the levels above 1% [1, 2]. Using SERS technique can significantly reduce the detection level for biopolymers [3, 4, 6]. By using a novel system for HTS Raman analysis, Identity plate reader, we studied two structurally different proteins: BSA, a globular protein with main structure represented by  $\alpha$ -helix and gelatin, a randomly coiled protein, by both conventional Raman spectroscopy (Figure 2) and SERS technique (Figure 3). Several anions at different concentrations were used as coagulation agents.


## APPLICATION NOTE

The results show that:

- 1). Identity plate reader can detect as low as 5% gelatin (Figure 2), and 1% BSA (not shown) in aqueous solutions by using conventional Raman technique;
- 2). using SERS technique significantly enhance sensitivity of protein detection (BSA at levels below 0.1% is easily identifiable);
- 3). chloride at final concentration above of 50-100mM was sufficient for the observed enhancement effect. Using phosphate, or sulfate at 100mM was as effective as that of chloride, while carbonate anion at neutral pH was totally ineffective;
- 4). SERS signal linearly depends on the protein concentration in a buffer containing 100mM chloride when the protein is present at concentrations below 1%.



**Figure 2.** Raman spectra of some proteins in water assayed by Identity equipped with 785nm Laser.



**Figure 3.** SERS (solid lines) and conventional Raman (dashed lines) spectra of BSA at different concentrations (0.2-5%) in 0.1M NaCl, with 1/10 diluted silver colloid. (Background subtracted; Y values paired up.)

## APPLICATION NOTE

**Recommended Procedure:** Switch Identity Raman plate reader on and start the Identity controlling software on the PC. Choose the “User Mode” for operation. For a system equipped with 785nm laser, we recommend using 100% power, while for the system equipped with 532nm laser, a 50% power is recommended for Raman data collection. Check inside the box for dark current correction. For detection of a protein in solution in the range of concentrations below 1%, add concentrated NaCl to the sample to make the final concentration of NaCl 100mM, dilute silver colloid with 100mM NaCl buffer 5-20 fold (the exact dilution factor needs to be determined experimentally for each batch of the colloid) and mix the sample with the diluted colloid in the well of the microplate immediately before the measurement (for a 96-well microtiter plate, use 0.1-0.15ml of sample and colloid solutions per well; for a 384-well plate, use 0.05ml of each solution per well). Perform Raman spectra analysis on Identity Raman plate reader using desired exposure time (from 0.01 to 60 sec per a single exposure) and necessary number of repeated exposures to be used for averaging of the final spectra. Collected spectral characteristics will be saved in files marked with the date (in the format year-month-day followed by the time (24hr based), the experiment has been started) and ready to be used by a third-party application software.

**References:**

1. J. L. Lippert, D. Tyminski, P. J. Desmeules “Determination of the secondary structure of proteins by laser Raman spectroscopy” *J. Am. Chem. Soc.*, 98: 7075–7080 (1976).
2. Peticolas WL “Raman spectroscopy of DNA and proteins” *Methods Enzymol.*, 246: 389-416 (1995).
3. X. X. Hun, G.G. Huang, B. Zhao and Y. Ozaki “Label-Free Highly Sensitive Detection of Proteins in Aqueous Solutions Using Surface-Enhanced Raman Scattering” *Anal. Chem.*, 81: 3329-3333 (2009).
4. Igor Zlatkin, Sarah Mendelowitz, Stuart Farquharson, David Drapcho. “Analysis of DNA-Polyamine Interactions Using SERS Method in Identity Raman Microplate Reader” *XXII International Conference on Raman Spectroscopy*, Boston, MA; Abstracts, Poster #099 (2010).
5. P. C. Lee and D. Meisel, “Adsorption and Surface-Enhanced Raman of Dyes on Silver and Gold Sols” *J. Phys. Chem.*, 86: 3391-3395 (1982).
6. D. L. Drapcho, I. Zlatkin, F. Inscore, C. Shende, A. Sengupta, H. Huang, and S. Farquharson, “High-Throughput Trace Analysis Using SERS-Active Microtiter Plates with a Raman Plate Reader” *Spectroscopy*, 25 (Suppl.): 42-50 (2010).

Please contact Dr. Igor Zlatkin ([izlatkin@digilabglobal.com](mailto:izlatkin@digilabglobal.com)) for additional technical information.

**WORLDWIDE OFFICE**  
84 October Hill Road  
Holliston, MA 01746  
United States  
Phone: 508 893 3130  
Fax: 508 893 8011

**EUROPEAN OFFICE**  
1B Blackstone Road  
Huntingdon, Cambridgeshire  
PE29 6EF United Kingdom  
Phone: [+44] 1480 426 700  
Fax: [+44] 1480 426 767

**ASIAN PACIFIC OFFICE**  
6th Fl. Yokohama World Porters  
2-2-1 Shinkou, Nakaku  
Yokohama, Japan 231-0001  
Phone/Fax: 045 651 6252